

SOMATOM Force - Scan protocols *syngo* CT VB20

The dose values stated in the scan protocol list are estimations and provide orientation and overview. For individual scan protocols, the valid dose values are displayed at the CT scanner.

The descriptions below refer to the corresponding column headings of the scan protocol tables on the following pages.

Qual. ref. tube current [Qual. ref. mAs]

'Qual. ref. mAs/rot.' for cardio scans

'Quality ref. mAs' indicates the preset reference values for operation with activated CARE Dose4D.

Tube current A [mAs]

'mAs' for sequential scans, 'eff. mAs' for spiral scans, 'mAs/rot.' for cardio scans

'mAs' indicates the preset reference values for operation with deactivated CARE Dose4D.

CTDIvol [mGy]

'mGy/s' for i-Fluoro scans

Measured in a 16 cm phantom for head modes and in a 32 cm phantom for body scans.

CTDIvol [mGy]; DLP [mGy cm];

The indicated dose values 'CTDIvol' and 'DLP' refer to the indicated settings (tube current A) for operation with deactivated CARE Dose4D, regardless if CARE Dose 4D is activated by default or not.

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
Head_XCARE	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Head		on	on	120		120None/None	332	332		47,32	649,87	16	standard		5,8	1	192	96	0,6	5	0,55		120	5	Hr40	
HeadRoutine	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Head		on	on	120		120None/None	332	332		47,32	649,87	16	standard		5,8	1	192	96	0,6	5	0,55		120	5	Hr40	
HeadNeuro_XCARE	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Head		on	on	120		120None/None	273	273		41,36	559,53	16	standard		7,69	1	128	64	0,6	5	0,55		120	5	Hr40	
HeadNeuro	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Head		on	on	120		120None/None	273	273		41,36	559,53	16	standard		7,69	1	128	64	0,6	5	0,55		120	5	Hr40	
HeadRoutineSeq	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	HeadSeq		on	on	120		120None/None	273	273		46,69	672,34	16	standard		1	1	192	96	0,6	5		48	139,3	5	Hr40	
HeadNeuroSeq	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	HeadSeq		on	on	120		120None/None	273	273		46,71	635,31	16	standard		1	1	128	64	0,6	5		34	133	5	Hr40	
InnerEarUHR	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr60	
	InnerEarUHR		on	on	130		130None/None	114	114		18,93	110,12	16	narrow		4,45	1	64	32	0,6	2	0,85		40	2	Ur77	
																					0,4				0,3	Ur77	
																					0,4				0,3	Ur77	
InnerEarUHR_Sn	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr60	
	InnerEarUHR		on	on	150		150TinFilter/None	473	473		25,79	121,51	16	narrow		10,33	1	64	32	0,6	2	0,25		40	2	Ur77	
																					0,4				0,3	Ur77	
																					0,4				0,3	Ur77	
InnerEar	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr60	
	InnerEar		on	on	130		130None/None	114	114		19,56	173,36	16	standard		3,26	1	192	96	0,6	2	0,8		58	2	Hr59	
																					0,6				0,4	Hr59	
																					0,6				0,4	Hr59	
InnerEarSeq	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr60	
	InnEarSeq		on	on	130		130None/None	114	114		19,56	112,64	16	standard		1	1	192	96	0,6	2		48	46,3	2	Hr59	
																					0,6				46,3	0,4	Hr59
Sinus	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Sinus		on	on	120		120None/None	35	35		4,99	55,19	16	standard		3,74	1	192	96	0,6	5	0,8		80	5	Hr59	
																					5			80	5	Hr40	
																					1			80	0,7	Hr59	
																					1			80	0,7	Hr40	
Sinus_Sn	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Sinus		on	on	100		100TinFilter/None	533	533		4,44	49,16	16	standard		3,74	1	192	96	0,6	5	0,8		80	5	Hr59	
																					5			80	5	Hr40	
																					1			80	0,7	Hr59	
																					1			80	0,7	Hr40	
Sinus_Vol	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Sinus		on	on	120		120None/None	35	35		4,99	55,19	16	standard		3,74	1	192	96	0,6	5	0,8		80	5	Hr59	
																					3				3	Hr40	
																					3				3	Hr59	
Orbit	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Orbit		on	on	120		120None/None	88	88		12,54	111,18	16	standard		3,26	1	192	96	0,6	5	0,8		58	5	Hr59	
																					5			58	5	Hr40	
																					0,6			58	0,4	Hr59	
Dental	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Dental		on	on	120		120None/None	63	63		8,98	79,6	16	standard		3,26	1	192	96	0,6	0,75	0,8		58	0,5	Hr59	
NeckRoutine	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20	
	Neck		on	on	120		120None/None	116	116		7,75	163,36	32	standard		5,91	1	192	96	0,6	5	0,8		180	5	Br40	
																					1			180	0,7	Br40	
NeckThorax	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Neck		on	on	120		120None/None	116	116		7,75	132,34	32	standard		5,04	1	192	96	0,6	5	0,8		140	5	Br40	
																					1			140	0,7	Br40	
		Pause																									
	Thorax		on	on	110		110None/None	51	51		2,66	88,02	32	standard		4,26	0,5	192	96	0,6	5	0,8		300	5	Br40	

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
TurboFlash_Thorax_HR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	FI_Thorax		on	on	110	110	110None/110None	110	110		4,61	168,1	32	narrow	narrow	0,93	0,25	96	96	0,6	5	3		300	5	Br40	
																					5			300	5	BI57	
																					0,75			300	0,5	Br40	
Thorax_XCARE	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Thorax		on	on	110		110None/None	51	51		2,66	85,19	32	standard		2,67	0,25	192	96	0,6	5	0,6		300	5	Br40	
																					5			300	5	BI57	
																					0,6			300	0,4	Br40	
ThoraxRoutine	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Thorax		on	on	110		110None/None	51	51		2,66	93,31	32	standard		3,17	0,5	192	96	0,6	5	1,2		300	5	Br40	
																					5			300	5	BI57	
																					0,6			300	0,4	Br40	
DSXXL_Thorax	Topogram				150		150TinFilter/None		20		0,02	0,83	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_Thorax		on	on	140	140	140None/140None	36	36		3,51	121,69	32	standard	standard	4,26	0,5	192	96	0,6	5	0,8		300	5	Br40	
																					5			300	5	BI57	
																					1			300	0,7	Br40	
Thorax_HR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	ThorHR		on	on	110		110None/None	110	110		5,63	197,5	32	standard		1,59	0,25	96	96	0,6	5	1,2		300	5	Br40	
																					5			300	5	BI57	
																					0,75			300	0,5	Br40	
ThoraxSeqHR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	ThorHRSeq		on	on	120		120None/None	110	110		1,36	40,94	32	standard		0,5	0,5	2	2	1	1		10	291		BI57	
	TurboFlash_LungLowDose_HR_Sn	Topogram			100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	FI_LungLowDo.		on	on	100	100	100TinFilter/100TinFilter	106	106		0,37	13,38	32	narrow	narrow	0,93	0,25	96	96	0,6	5	3		300	5	Br40	
																					5			300	5	BI57	
																					1			300	0,7	Br40	
																					1			300	0,7	Br64	
																								512		Tr20	
LungLowDose_Sn	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	LungLowDose		on	on	100		100TinFilter/None	101	101		0,43	15	32	standard		3,17	0,5	192	96	0,6	5	1,2		300	5	Br40	
																					5			300	5	BI57	
																					1			300	0,7	Br40	
LungLowDose_HR_Sn	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	LungLowDose		on	on	100		100TinFilter/None	101	101		0,43	15	32	standard		1,59	0,25	96	96	0,6	5	1,2		300	5	Br40	
																					5			300	5	BI57	
																					1			300	0,7	Br40	
LungCAD	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	LungCAD		on	on	120		120None/None	20	20		1,34	46,87	32	standard		3,17	0,5	192	96	0,6	5	1,2		300	5	Br40	
																					5			300	5	BI57	
																					1			300	0,7	Br40	
ThorAbd_XCARE	Topogram				100		100TinFilter/None		75		0,01	0,68	32	standard		7,85		6	6	0,6				768		Tr20	
	ThorAbd		on	on	120		120None/None	147	147		9,83	609,34	32	standard		9,7	0,5	192	96	0,6	5	0,6		600	5	Br40	
																					5				5	BI57	
																					5				5	Br40	
ThorAbd	Topogram				100		100TinFilter/None		75		0,01	0,68	32	standard		7,85		6	6	0,6				768		Tr20	
	ThorAbd		on	on	120		120None/None	147	147		9,83	609,34	32	standard		9,7	0,5	192	96	0,6	5	0,6		600	5	Br40	
																					5				5	BI57	
																					5				5	Br40	
ThoraxIntervention	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	i-Spiral		on	on	110		110None/None	138	138		7,21	252,68	32	standard		3,17	0,5	48	48	1,2	5	1,2		300	3,5	Br40	
	i-Fluoro		off	off	110		110None/None	38	38		7,02	10,11	32	standard		99	0,5	12	12	1,2	4,8		0	9,6	4,8	Br40	
	i-Sequence		off	off	110		110None/None	38	38		2,54	3,65	32	standard		0,5	0,5	12	12	1,2	4,8		0	9,6		Br40	
AbdomenRoutine	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Abdomen		on	on	120		120None/None	147	147		9,83	216,25	32	standard		3,9	0,5	192	96	0,6	5	0,6		200	5	Br40	
																					1			200	0,7	Br40	

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
AbdMultiPhase	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Non Contrast		on	on	120		120None/None	133	133		8,9	195,84	32	standard		3,9	0,5	48	48	1,2	5	0,6		200	5	Br40	
		Pause																									
	Arterial Phase		on	on	120		120None/None	147	147		9,83	216,25	32	standard		3,9	0,5	192	96	0,6	5	0,6		200	5	Br40	
																					1			200	0,7	Br40	
AbdomenSeq	Venous Phase		on	on	120		120None/None	133	133		8,9	195,84	32	standard		3,9	0,5	48	48	1,2	5	0,6		200	5	Br40	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	AbdSeq		on	on	120		120None/None	147	147		11,79	283,03	32	standard		0,5	0,5	192	96	0,6	5		48	235,3	5	Br40	
ColonLowDose_Sn	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Colo_supine		on	on	150		150TinFilter/None	79	79		2,3	100,37	32	standard		4,86	0,5	192	96	0,6	1	0,9		400	0,7	Br40	
		Pause PR																									
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
ColonCAD	Colo_prone		on	on	150		150TinFilter/None	43	43		1,25	56,67	32	standard		3,48	0,5	192	96	0,6	1	1,4		400	0,7	Br36	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Colo_supine		on	semi	120		120None/None	55	55		3,68	160,25	32	standard		4,86	0,5	192	96	0,6	1	0,9		400	0,7	Br40	
		Pause PR																									
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
WholeBodyBone_Sn	Colo_prone		on	semi	120		120None/None	30	30		2,01	90,68	32	standard		3,48	0,5	192	96	0,6	1	1,4		400	0,7	Br36	
	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				1024		Tr20	
	WholeBody		on	on	100		100TinFilter/None	653	653		2,77	226,92	32	standard		25,19	1	192	96	0,6	5	0,6		800	5	Br59	
AbdomenIntervention																					1			800	0,7	Br59	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	i-Spiral		on	on	120		120None/None	210	210		14,05	464,6	32	standard		4,26	0,5	48	48	1,2	5	0,8		300	3,5	Br40	
PelvisRoutine	i-Fluoro		off	off	120		120None/None	50	50		11,84	17,04	32	standard		99	0,5	12	12	1,2	4,8		0	9,6	4,8	Br40	
	i-Sequence		off	off	120		120None/None	50	50		4,27	6,15	32	standard		0,5	0,5	12	12	1,2	4,8		0	9,6		Br40	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr60	
DSXXL_Pelvis	Pelvis		on	on	120		120None/None	147	147		9,84	226,87	32	standard		3,17	0,5	48	48	1,2	5	0,8		200	5	Br40	
	Topogram				150		150TinFilter/None		20		0,02	0,83	32	standard		5,29		6	6	0,6				512		Tr60	
	DS Pelvis		on	on	140	140	140None/140None	120	120		11,7	269,86	32	standard	standard	3,9	0,5	192	96	0,6	5	0,6		200	5	Br40	
Hip	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr60	
	Hip		on	on	120		120None/None	133	133		8,89	116,16	32	standard		4,17	1	192	96	0,6	5	0,8		100	5	Br59	
																					0,6			100	0,4	Br59	
Hip_iMAR	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr60	
	Hip		on	on	120		120None/None	133	133		8,89	116,16	32	standard		4,17	1	192	96	0,6	5	0,8		100	5	Br59	
																								100	5	Br59	
																					0,6			100	0,4	Br59	
																					0,6			100	0,4	Br59	
C_Spine_Vol	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr60	
	C_Spine		on	on	120		120None/None	193	193		12,9	245,99	32	standard		5,48	1	192	96	0,6	2	0,8		160	2	Br40	
																					0,6				0,4	Br40	
																					2				2	Br40	
																					2				2	Br59	
C_Spine	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr60	
	C_Spine		on	on	120		120None/None	193	193		12,9	245,99	32	standard		5,48	1	192	96	0,6	2	0,8		160	2	Br40	
																					0,6			160	0,4	Br40	
																					0,6			16			

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
SpineNeuro Vol	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr60	
	Spine		on	on	130		130None/None	193	193		16,75	308,72	32	standard		7,21	1	128	64	0,6	3	0,8		160	3	Br40	
																				0,6					0,4	Br40	
																				2					2	Br40	
																				2					2	Br59	
SpineNeuro	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr60	
	Spine		on	on	130		130None/None	193	193		16,75	308,72	32	standard		7,21	1	128	64	0,6	3	0,8		160	3	Br40	
																				0,6				160	0,4	Br40	
																				0,6				160	0,4	Br59	
SpineSeq	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr60	
	L3-4		on	on	130		130None/None	129	129		13	24,96	32	standard		1	1	64	32	0,6	2		17,5	15,2	2	Br40	
	L4-5		on	on	130		130None/None	129	129		13	24,96	32	standard		1	1	64	32	0,6	2		17,5	15,2	2	Br40	
	L5-S1		on	on	130		130None/None	129	129		13	24,96	32	standard		1	1	64	32	0,6	2		17,5	15,2	2	Br40	
SpineIntervention	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20	
	i-Spiral		on	on	130		130None/None	230	230		18,8	358,39	32	standard		2,74	0,5	48	48	1,2	5	0,8		160	3,5	Br40	
	i-Fluoro		off	off	120		120None/None	60	60		14,2	20,45	32	standard		99	0,5	12	12	1,2	4,8		0	9,6	4,8	Br40	
	i-Sequence		off	off	120		120None/None	60	60		5,13	7,39	32	standard		0,5	0,5	12	12	1,2	4,8		0	9,6		Br40	
SpineBoneIntervention_Sn	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20	
	i-Spiral		on	on	100		100TinFilter/None	1598	1598		6,78	112,76	32	standard		3,78	0,5	48	48	1,2	5	0,35		160	3,5	Br59	
	i-Fluoro		off	off	100		100TinFilter/None	266	266		3,99	5,75	32	standard		99	0,5	12	12	1,2	4,8		0	9,6	4,8	Br64	
	i-Sequence		off	off	100		100TinFilter/None	400	400		2,17	3,12	32	standard		0,5	0,5	12	12	1,2	4,8		0	9,6		Br64	
WristUHR	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	WristUHR		on	off	120		120None/None	84	84		4,87	38,06	32	narrow		5,68	1	64	32	0,6	2	0,85		60	2	Ur77	
																				0,4				60	0,2	Ur77	
ExtrUHR	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	ExtrUHR		on	off	120		120None/None	84	84		4,87	38,06	32	narrow		5,68	1	64	32	0,6	2	0,85		60	2	Ur77	
																				0,4				60	0,3	Ur77	
ExtrUHR_Sn	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	ExtrUHR		on	off	150		150TinFilter/None	220	220		5,82	42,35	32	narrow		7,71	1	64	32	0,6	2	0,55		60	2	Ur77	
																				0,4				60	0,3	Ur77	
Extremity	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	Extremity		on	off	120		120None/None	70	70		4,68	84,54	32	standard		5,26	1	192	96	0,6	4	0,8		150	4	Br59	
																				1				150	0,7	Br59	
Extremity_Sn	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	Extremity		on	off	100		100TinFilter/None	932	932		3,95	71,35	32	standard		5,26	1	192	96	0,6	4	0,8		150	4	Br59	
																				1				150	0,7	Br59	
KneeUHR	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	KneeUHR		on	off	120		120None/None	98	98		5,68	95,51	32	narrow		11,2	1	64	32	0,6	3	0,85		150	3	Ur77	
																				0,4				150	0,2	Ur77	
FootUHR	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20	
	FootUHR		on	off	120		120None/None	84	84		4,87	81,87	32	narrow		11,2	1	64	32	0,6	3	0,85		150	3	Ur77	
																				0,4				150	0,2	Ur77	
ExtrUHR	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	ExtrUHR		on	off	120		120None/None	84	84		4,87	81,87	32	narrow		11,2	1	64	32	0,6	2	0,85		150	2	Ur77	
																				0,4				150	0,3	Ur77	
ExtrUHR_Sn	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	ExtrUHR		on	off	150		150TinFilter/None	220	220		5,82	94,73	32	narrow		16,29	1	64	32	0,6	2	0,55		150	2	Ur77	
																				0,4				150	0,3	Ur77	
Extremity	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	Extremity		on	off	120		120None/None	70	70		4,68	84,54	32	standard		5,26	1	192	96	0,6	4	0,8		150	4	Br59	
																				1				150	0,7	Br59	
Extremity_Sn	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	Extremity		on	off	100		100TinFilter/None	932	932		3,95	71,35	32	standard		5,26	1	192	96	0,6	4	0,8		150	4	Br59	
																				1				150	0,7	Br59	
HeadAngio_XCARE	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	HeadAngio		on	on	90		90None/None	133	133		9,26	126,94	16	standard		1,8	0,25	128	64	0,6	4	0,6		120	4	Hv36	
																				0,6				120	0,4	Hv36	
HeadAngioRoutine	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	HeadAngio		on	on	90		90None/None	133	133		9,26	145,95	16	standard		2,3	0,5	128	64	0,6	4	1,2		120	4	Hv36	
																				0,6				120	0,4	Hv36	
HeadAngio_Vol	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	HeadAngio		on	on	90		90None/None	133	133		9,26	145,95	16	standard		2,3	0,5	128	64	0,6	4	1,2		120	4	Hv36	
																				4					4	Hv36	
HeadAngioUltraFastMultiple	Topogram				150		150TinFilter/None		50		0,08	3,85	16	standard		5,29		6	6	0,6				512		Tr20	
	HeadAngioUFM		off	off	90		90None/None	133	133		18,93	810,86	16	standard		2,47	0,25	192	96	0,6	4	1	48	454	4	Hv36	
																				0,6				454	0,4	Hv36	
NeuroDSACT	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	Non Contrast		on	on	90		90None/None	86	86		5,99	94,37	16	standard		2,3	0,5	128	64	0,6	0,6	1,2		120	0,4	Hv36	
		Pause																									
	HeadAngio		on	on	90		90None/None	133	133		9,26	145,95	16	standard		2,3	0,5	128	64	0,6	4	1,2		120	4	Hv36	
																				0,6				120	0,4	Hv36	
CarotidAngio	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	CarotidAngio		on	on	120																						

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
CarotidAngio_Vol	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	CarotidAngio		on	on	120		120None/None	84	84		5,97	201,56	32	standard		2,13	0,25	128	64	0,6	5	1,2		300	5	Bv36	
																					3				3	Bv36	
CarotidDSACT	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	Non Contrast		on	on	120		120None/None	42	42		2,98	100,78	32	standard		2,13	0,25	128	64	0,6	0,6	1,2		300	0,4	Bv40	
		Pause																									
	CarotidAngio		on	on	120		120None/None	84	84		5,97	201,56	32	standard		2,13	0,25	128	64	0,6	5	1,2		300	5	Bv36	
																					0,6			300	0,4	Bv40	
TurboFlash_ThorAngio	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	FI_ThorAngio		on	on	80	80	80None/80None	174	174		2,55	116,49	32	narrow	narrow	1,41	0,25	192	96	0,6	5	1,9		400	5	Bv36	
																					5			400	5	Bv59	
																					0,6			400	0,4	Bv36	
ThorAngio_XCARE	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	ThorAngio		on	on	80		80None/None	157	157		3,01	126,63	32	standard		3,4	0,25	192	96	0,6	5	0,6		400	5	Bv36	
																					5			400	5	Bv59	
																					0,6			400	0,4	Bv36	
ThorAngioRoutine	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	ThorAngio		on	on	80		80None/None	157	157		3,01	135,81	32	standard		3,89	0,5	192	96	0,6	5	1,2		400	5	Bv36	
																					5			400	5	Bv59	
																					0,6			400	0,4	Bv36	
ThorAngio_Vol	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	ThorAngio		on	on	80		80None/None	157	157		3,01	135,81	32	standard		3,89	0,5	192	96	0,6	5	1,2		400	5	Bv36	
																					5			400	5	Bv59	
																					3				3	Bv36	
TurboFlash_ChestPainECG	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
		Contrast																									
	TestBolus		off	off	100		100None/None	24	24		12,94	12,94	32	standard		0,25	0,25	1	1	10	10		0	0		Bv40	
		Contrast																									
	FI_ChestPain		on	on	100	100	100None/100None	288	288		2,85	83,09	32	narrow	narrow	0,84	0,25	192	96	0,6	3	3,2		250	3	Bv36	
																					0,6				0,3	Bv36	
																					3			250	3	Bv36	
																					3				3	Bv59	
																					1				0,5	Bv40	
		Pause																									
	DS_CorAdSeq		on	on	100	100	100None/100None	288	288		4,96	71,49	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48	150,6	3	Bv36	
																					0,6			150,6	0,4	Bv36	
DS_ChestPainECG	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_ChestPain		on	on	100	100	100None/100None	288	288		8,18	224,94	32	narrow	narrow	6,94	0,25	192	96	0,6	3	0,17		250	3	Bv36	
																					0,6				0,3	Bv36	
																					0,6				0,3	Bv36	
																					3			250	3	Bv36	
																					3				3	Bv59	
																					3				3	Bv36	
DS_ChestECG_Abd	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				1024		Tr20	
	DS_ChestPain		on	on	100	100	100None/100None	288	288		8,18	224,94	32	narrow	narrow	6,94	0,25	192	96	0,6	3	0,17		250	3	Bv36	
																					0,6				0,3	Bv36	
																					0,6				0,3	Bv36	
																					3			250	3	Bv36	
																					3				3	Bv59	
																					3				3	Bv36	
																					1				0,5	Bv36	
DS_ChestECG_Abd	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				1024		Tr20	
	DS_ChestPain		on	on	100	100</																					

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
TurboFlash_CaScore	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	FI_CaSc		on	on	120	120	120None/120None	80	80		1,33	23,57	32	narrow	narrow	0,68	0,25	192	96	0,6	3	3,2		135	1,5	Sa36	
DS_CaScore	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CaSc		on	on	120	120	120None/120None	80	80		3,83	67,06	32	narrow	narrow	4,37	0,25	192	96	0,6	3	0,17		150	1,5	Sa36	
DS_CaScoreSeq_Sn	Topogram				100	100	100TinFilter/100None		75		0,01	0,45	32	standard	standard	5,29		6	6	0,6				512		Tr20	
DS_CaScoreSeq	DS_CaScSeq		on	on	100	100	100TinFilter/100TinFilter	534	534		1,6	23,07	32	narrow	narrow	0,18	0,25	48	48	1,2	3		48	150,6	1,5	Sa36	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
DS_CaScSeq	DS_CaScSeq		on	on	120	120	120None/120None	80	80		3,7	53,35	32	narrow	narrow	0,18	0,25	48	48	1,2	3		48	150,6	1,5	Sa36	
FAST_TurboFlash_Cardiac	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	FI_CaSc		on	on	120	120	120None/120None	80	80		1,33	23,57	32	narrow	narrow	0,68	0,25	192	96	0,6	3	3,2		135	1,5	Sa36	
	TestBolus	Contrast	off	off	100		100None/None	24	24		12,94	12,94	32	standard		0,25	0,25	1	1	10	10		0	0		Bv40	
	FI_CorCTA	Contrast	on	on	100	100	100None/100None	288	288		2,85	50,32	32	narrow	narrow	0,68	0,25	192	96	0,6	3	3,2		135	3	Bv40	
		Pause																			0,6			135	0,3	Bv40	
	DS_CorAdSeq		on	on	100	100	100None/100None	288	288		4,96	71,49	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48	150,6	3	Bv40	
																					0,6			150,6	0,4	Bv40	
DS_CoroCTA_AdaptSeq	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorAdSeq		on	on	100	100	100None/100None	288	288		4,96	71,49	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48	150,6	3	Bv40	
																					0,6			150,6	0,4	Bv40	
DS_CoroCTA_AdaptSeq_Edge_IR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorAdSeq		on	on	100	100	100None/100None	288	288		4,96	71,49	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48	150,6	3	Bv40	
																					0,6			150,6	0,4	Bv40	
																					0,5			150,6	0,15	Bv49	
DS_CoronaryCTA	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorCTA		on	on	100	100	100None/100None	288	288		17,49	253,88	32	narrow	narrow	3,59	0,25	192	96	0,6	3	0,17		120	3	Bv40	
																					0,6			120	0,3	Bv40	
																					1,5			120	1	Qr36	
DS_CoroCTA_BiSeg	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorCTABi		on	on	100	100	100None/100None	288	288		17,49	253,88	32	narrow	narrow	3,59	0,25	192	96	0,6	3	0,17		120	3	Bv40	
																					0,6			120	0,4	Bv40	
																					1,5			120	1	Qr36	
DS_CoronaryCTA05s	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorCTA05s		on	on	100	100	100None/100None	288	288		17,49	253,88	32	narrow	narrow	7,19	0,5	192	96	0,6	3	0,17		120	3	Bv40	
																					0,6			120	0,4	Bv40	
																					1,5			120	1	Qr36	
DS_CoroCTA_AdaptSeqArrhythmicHR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorAdSeq		on	on	100	100	100None/100None	288	288		13,13	189,12	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48	150,6	3	Bv40	
																					0,6			150,6	0,4	Bv40	
DS_CoroCTAArrhythmicHR	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorCTA		on	on	100	100	100None/100None	288	288		16,09	233,56	32	narrow	narrow	3,59	0,25	192	96	0,6	3	0,17		120	3	Bv40	
																					0,6			120	0,3	Bv40	
DS_CoroCTA_HeartFlow	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DS_CorAdSeq		on	on	100	100	100None/100None	340	340		5,86	84,39	32	narrow	narrow	0,18	0,25	192	96	0,6	0,6		48	153	0,4	Bv40	
																					0,6			153	0,4	Bv40	
TurboFlash_TAVI	Topogram				100	100	100TinFilter/100None		75		0,01	0,68	32	standard	standard	7,85		6	6								

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
DE AbdomenDynSerio	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	DE AbdSeqDyn		on	off	90	150	90None/150TinFilter	215	215	134	10,11	58,25	32	standard	standard	0,5	0,5	192	96	0,6	5		0	52,6	5	Br40	
DE Body BoneMarrow	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				52,6	1	Qr40	
	DE Body		on	off	90	150	90None/150TinFilter	220	220	138	11,03	798,78	32	standard	standard	16,22	0,5	128	64	0,6	5	0,6		700	5	Br40	
DE Spine Metal Vol	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				700	1	Qr40	
	DE Spine		on	off	100	150	100None/150TinFilter	190	190	380	19,87	356,78	32	standard	standard	5,17	0,5	128	64	0,6	3	0,5		512	3	Br40	
																					1,5			160	1	Qr40	
																					1,5			160	1	Br40	
																					2				2	2	Br40
																					2				2	2	Br59
DE Hand Gout	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	DE Hand		on	off	80	150	80None/150TinFilter	100	100	67	4,12	73,68	32	standard	standard	3,8	0,5	128	64	0,6	2	0,7		150	2	Br40	
DE Foot Gout	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				150	0,5	Qr40	
	DE Foot		on	off	80	150	80None/150TinFilter	125	125	83	5,12	91,68	32	standard	standard	3,8	0,5	128	64	0,6	2	0,7		256	2	Br40	
DE Knee BoneMarrow	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				150	0,5	Qr40	
	DE Knee		on	off	80	150	80None/150TinFilter	150	150	100	7,15	121,71	32	standard	standard	6,6	0,5	64	32	0,6	3	0,7		256	3	Br40	
DE HeadAngio BoneRem	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				150	0,5	Qr40	
	DE HeadAngio		on	off	80	150	80None/150TinFilter	100	100	67	10,02	140,41	16	standard	standard	2,74	0,25	64	32	0,6	4	0,7		256	4	Tr20	
DE Head BrainHem	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				120	0,7	Qr40	
	DE Head		on	off	80	150	80None/150TinFilter	310	310	207	31	434,63	16	standard	standard	10,96	1	64	32	0,6	4	0,7		120	4	Hr40	
DE CarotidAngio BoneRem	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				120	0,7	Qr40	
	DE CarotidAngio		on	off	90	150	90None/150TinFilter	90	90	69	3,95	110,16	32	narrow	narrow	2,83	0,25	128	64	0,6	5	0,7		250	5	Bv36	
DE CarotidAngioFast BoneRem	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				250	0,7	Qr40	
	DE CarotidAngio		on	off	90	150	90None/150TinFilter	90	90	69	4,61	132,55	32	standard	standard	2,05	0,25	192	96	0,6	5	0,7		256	5	Bv36	
DE BodyAngio BoneRem	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				250	0,7	Qr40	
	DE BodyAngio		on	off	90	150	90None/150TinFilter	95	95	59	4,74	203,49	32	standard	standard	8,46	0,5	128	64	0,6	5	0,7		512	5	Tr20	
DE BodyAngioFast BoneRem	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				400	1	Bv36	
	DE BodyAngio		on	off	90	150	90None/150TinFilter	95	95	59	4,46	195,26	32	standard	standard	2,98	0,25	192	96	0,6	5	0,7		400	1	Qr40	
DE ExtrAngio BoneRem	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				400	1	Qr40	
	DE ExtrAngio		on	off	80	150	80None/150TinFilter	120	120	67	4,53	465,66	32	standard	standard	9,83	0,25	128	64	0,6	5	0,7		1024	5	Bv36	
DE CoroCTA HeartPBV	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				1000	1	Qr40	
	DE CorCTA		on	off	90	150	90None/150TinFilter	165	165	127	19,48	275,72	32	narrow	narrow	5,12	0,25	128	64	0,6	3	0,17		512	3	Tr20	
DE CoroCTA AdaptSeq	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				120	0,3	Bv40	
	DE CorAdSeq		on	off	90	150	90None/150TinFilter	165	165	127	5,43	139,3	32	narrow	narrow	0,25	0,25	128	64	0,6	3		34	134,5	3	Qr40	
																					0,6			120	0,7	Tr20	
																					1			120	0,7	Tr20	
RDE Thorax Analysis XCARE	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				134,5	0,4	Bv40	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				134,5	1	Qr40	
DE Thorax	Topogram		on	off	90	150	90None/150TinFilter	110	110	85	5,65	184,96	32	standard	standard	2,87	0,25	192	96	0,6	5	0,55		512	5	Tr20	
																					5			300	5	Bl64	
RDE Thorax Analysis	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				300	1	Qr40	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
DE Thorax	Topogram		on	off	90	150	90None/150TinFilter	110	110	85	5,65	184,96	32	standard	standard	2,87	0,25	192	96	0,6	5	0,55		300	5	Br40	
																					5			300	5	Bl64	
RDE Abdomen_VNC	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				300	1	Qr40	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
DE Abdomen	Topogram		on	off	100	150	100None/150TinFilter	240	240	120	13,95	312,83	32	standard	standard	5,35	0,5	128	64	0,6	5	0,6		200	5	Br40	
																					1,5			200	1	Qr40	
RDE AbdomenSeq	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
DE AbdSeq	Topogram		on	off	100	150	100None/150TinFilter	190	190	95	12,47	254,42	32	standard	standard	0,5	0,5	128	64	0,6	5		34	202	5	Br40	
																					1,5			202	1	Qr40	
RDE AbdomenDynSerio	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20	
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
DE AbdSeqDyn	Topogram		on	off	90	150	90None/150TinFilter	215	215	134	10,11	58,25	32	standard	standard	0,5	0,5	192	96	0,6	5		0	52,6	5	Br40	
																					1,5			52,6	1	Qr40	

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel		
Adult																												
RDE Body BoneMarrow	Topogram				150		150TinFilter/None		50		0,04	4,15	32	standard		10,41		6	6	0,6				1024		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard		10,41		6	6	0,6				1024		Tr20		
	DE Body		on	off	90	150	90None/150TinFilter	280	280	175	14,02	1015,2	32	standard	standard	16,22	0,5	128	64	0,6	5	0,6		700	5	Br40		
																					1,5			700	1	Qr40		
RDE Spine Metal Vol	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20		
	DE Spine		on	off	100	150	100None/150TinFilter	190	190	380	19,87	356,78	32	standard	standard	5,17	0,5	128	64	0,6	3	0,5		160	3	Br40		
																					1,5			160	1	Qr40		
																										1,5		Br40
																					2				2		Br40	
																					2				2		Br59	
RDE Hand Gout	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20		
	DE Hand		on	off	80	150	80None/150TinFilter	205	205	137	8,43	150,85	32	standard	standard	3,8	0,5	128	64	0,6	2	0,7		150	2	Br40		
																					0,75			150	0,5	Qr40		
RDE Foot Gout	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20		
	DE Foot		on	off	80	150	80None/150TinFilter	250	250	167	10,28	183,92	32	standard	standard	3,8	0,5	128	64	0,6	2	0,7		150	2	Br40		
																					0,75			150	0,5	Qr40		
RDE Knee BoneMarrow	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20		
	DE Knee		on	off	80	150	80None/150TinFilter	240	240	160	11,44	194,73	32	standard	standard	6,6	0,5	64	32	0,6	3	0,7		150	3	Br40		
																					0,75			150	0,5	Qr40		
RDE HeadAngio BoneRem	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,02	0,45	16	standard		2,73		6	6	0,6				256		Tr20		
	DE HeadAngio		on	off	80	150	80None/150TinFilter	155	155	103	15,47	216,85	16	standard	standard	2,74	0,25	64	32	0,6	4	0,7		120	4	Hv36		
																					1			120	0,7	Qr40		
RDE Head BrainHem	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,02	0,45	16	standard		2,73		6	6	0,6				256		Tr20		
	DE Head		on	off	80	150	80None/150TinFilter	310	310	207	31	434,63	16	standard	standard	10,96	1	64	32	0,6	4	0,7		120	4	Hr40		
																					1			120	0,7	Qr40		
RDE CarotidAngio BoneRem XCARE	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20		
	DE CarotidAngio		on	off	90	150	90None/150TinFilter	145	145	112	6,38	175,1	32	narrow	narrow	3,22	0,25	128	64	0,6	5	0,6		250	5	Bv36		
																					1			250	0,7	Qr40		
RDE CarotidAngio BoneRem	Topogram				150		150TinFilter/None		50		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20		
	DE CarotidAngio		on	off	90	150	90None/150TinFilter	145	145	112	6,38	178,08	32	narrow	narrow	2,83	0,25	128	64	0,6	5	0,7		250	5	Bv36		
																					1			250	0,7	Qr40		
RDE BodyAngio BoneRem	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20		
	DE BodyAngio		on	off	90	150	90None/150TinFilter	150	150	94	7,52	322,43	32	standard	standard	8,46	0,5	128	64	0,6	5	0,7		400	5	Bv36		
																					1,5			400	1	Qr40		
RDE Body OptiCon	Topogram				150		150TinFilter/None		50		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20		
	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20		
	DE Body		on	off	80	140	80None/140None	176	176	32	6,91	296,37	32	standard	standard	8,46	0,5	128	64	0,6	5	0,7		400	5	Bv36		

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
Trauma	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard			10,41	6	6	0,6					1024		Tr20
	Trauma		on	on	120		120None/None	180	180		12,04	966,97	32	standard			5,34	0,5	48	48	1,2	7	1,5		750	7	Br40
																						1,5			750	1	Br40
																						1,5			750	1	Br64
TraumaThin	Topogram				100		100TinFilter/None		75		0,01	0,91	32	standard			10,41	6	6	0,6					1024		Tr20
	TraumaThin		on	on	120		120None/None	180	180		12,79	1009,6	32	standard			7,51	0,5	128	64	0,6	7	1,5		750	7	Br40
																						0,6			750	0,3	Br40
																						0,6			750	0,3	Br64
PolyTrauma	Topogram				150		150TinFilter/None		50		0,08	3,85	16	standard			5,29	6	6	0,6					512		Tr20
	Head		on	on	120		120None/None	340	340		48,5	666,24	16	standard			2,9	0,5	48	48	1,2	6	0,55		120	6	Hr40
	Neck		on	on	120		120None/None	165	165		11,04	188,41	32	standard			2,52	0,5	48	48	1,2	5	0,8		140	5	Br40
		Pause																									
	Topogram				100		100TinFilter/None		75		0,01	0,68	32	standard			7,85	6	6	0,6					768		Tr20
	ThorAbd		on	on	120		120None/None	210	210		14,05	916,26	32	standard			4,72	0,5	48	48	1,2	7	1,4		600	7	Br40
																									7		BI57
																									7		Br40
																									1,5		Br40
TurboFlash_HeadTrauma	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard			2,73	6	6	0,6					256		Tr20
	FI_Head		on	on	120	120	120None/120None	210	210		29,93	518,72	16	standard	standard		0,84	0,25	192	96	0,6	5	1,55		120	5	Hr40
																									5		Hr59
ThoraxUHR_Sn	Topogram				100		100TinFilter/None		75		0,01	0,45	32	standard			5,29	6	6	0,6					512		Tr20
	ThorUHR		on	on	150		150TinFilter/None	91	91		2,24	71,86	32	narrow			8,81	0,5	32	32	0,6	5	1		300	5	UI64
																									5		Ub36
																									0,6		UI57
																									0,6		Ub36
NeuroVPCT_Stroke	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard			2,73	6	6	0,6					256		Tr20
	Head		on	on	120		120None/None	320	320		48,48	655,86	16	standard			7,69	1	128	64	0,6	5	0,55		120	5	Hr40
																									1		Hr36
		Pause																									
	DynMulti4D		off	off	70		70None/None	200	200		147,7	2183,19	16	standard			1,5	0,25	48	48	1,2	5	0,5	48	114	3	Hr36
																									1,5		Hr36
		Pause																									
	HeadAngio		on	on	90		90None/None	190	190		13,23	208,49	16	standard			2,3	0,5	128	64	0,6	1	1,2		120	0,7	Hr36
NeuroVPCT_Prolonge	Topogram				150		150TinFilter/None		50		0,08	1,92	16	standard			2,73	6	6	0,6					256		Tr20
	Head		on	on	120		120None/None	320	320		48,48	655,86	16	standard			7,69	1	128	64	0,6	5	0,55		120	5	Hr40
																									1		Hr36
		Pause																									
	DynMulti4D		off	off	70		70None/None	200	200		147,7	2183,19	16	standard			1,5	0,25	48	48	1,2	5	0,5	48	114	3	Hr36
																									1,5		Hr36
		Pause																									
	HeadAngio		on	on	90		90None/None	190	190		13,23	208,49	16	standard			2,3	0,5	128	64	0,6	1	1,2		120	0,7	Hr36
NeuroPCT	DynMulti		off	off	70		70None/None	200	200		229,19	1320,15	16	standard			40	1	48	48	1,2	5		48	52,6	3	Hr36
NeuroPCTSeq	DynSerio		off	off	70		70None/None	200	200		171,89	990,11	16	standard			0,5	0,5	48	48	1,2	5		0	52,6	3	Hr36
NeuroMultiPhase_VAngio	DynMulti4D		off	off	70		70None/None	200	200		177,4	3349,62	16	standard			1,5	0,325	48	48	1,2	5	1	48	176	3	Hr36
																									1,5		Hr36
NeuroMultiPhase_VAngio_long	DynMulti4D_L		off	off	80		80None/None	150	150		69,19	1913,35	16	standard			2	0,325	48	48	1,2	5	1	48	265	3	Hr36
																									1,5		Hr36
BodyVPCT	DynMulti4D		off	off	80		80None/None	150	150		74,31	1098,3	32	standard			1,5	0,25	48	48	1,2	5	0,5	48	114	3	Br36
																									3		Br36
																									1,5		Br36
		adjust to spec																									
BodyPCT	DynMulti		off	off	80		80None/None	150	150		230,6	1328,26	32	standard			40	0,5	48	48	1,2	5		48	52,6	3	Br36
		adjust to spec																									
BodyPCTSeq	DynSerio		off	off	80		80None/None	150	150		86,48	498,1	32	standard			0,5	0,5	48	48	1,2	5		0	52,6	3	Br36
		adjust to spec																									
BodyMultiPhase_VAngio	DynMulti4D_L		off	off	80		80None/None	150	150		29,52	816,37	32	standard			2	0,325	48	48	1,2	5	1	48	265	3	Br36
																									1,5		Br36
		adjust to spec																									
ExtrMultiPhase_VAngio	DynMulti4D_XL		off	off	80		80None/None	100	100		19,37	1150,63	32	standard			2,5	0,25	48	48	1,2	5	1,46	48	630	3	Br36
																									1,5		Br36
		adjust to spec																									
HeartVPCT	DynSerio4D		on	on	80	80	80None/80None	300	300		30,19	318,77	32	narrow	narrow		0,61	0,25	48	48	1,2	3		48	102,6	2	Qr36

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Adult																											
Sn_Topos	TopoHead				150		150TinFilter/None		137		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	TopoHeadNeck				150		150TinFilter/None		137		0,08	1,92	16	standard		2,73		6	6	0,6				256		Tr20	
	TopoNeck				150		150TinFilter/None		137		0,04	1,04	32	standard		2,73		6	6	0,6				256		Tr20	
	TopoNeckThor				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoShoulder				100		100TinFilter/None		205		0,01	0,23	32	standard		2,73		6	6	0,6				256		Tr20	
	TopoThorax				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoThorUpAbd				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoTorAbd				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoAbdomen				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoPelvis				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoLegs				100		100TinFilter/None		397		0,01	0,45	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoSpine				150		150TinFilter/None		265		0,04	2,08	32	standard		5,29		6	6	0,6				512		Tr20	
	TopoRunoff				100		100TinFilter/None		781		0,01	0,91	32	standard		10,41		6	6	0,6				1024		Tr20	
	TopoWholeBody				100		100TinFilter/None		1289		0,01	1,5	32	standard		17,19		6	6	0,6				1700		Tr20	

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Child																											
Head_XCARE	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Head		on	on	100		100None/None	495	248		22	308,14	16	standard		5,48	1	192	96	0,6	5	0,6		120	5	Hr40	
HeadRoutine	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Head		on	on	100		100None/None	495	248		22	331,47	16	standard		4,61	1	192	96	0,6	5	0,8		120	5	Hr40	
HeadNeuro	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Head		on	on	100		100None/None	407	204		19,24	277,59	16	standard		5,91	1	128	64	0,6	5	0,8		120	5	Hr40	
HeadRoutineSeq	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	HeadSeq		on	on	100		100None/None	407	204		21,72	417	16	standard		1	1	192	96	0,6	5		48	187,3	5	Hr40	
NeonateHead	Topogram				100		100TinFilter/None		20				16	narrow		2,73		6	6	0,6				256		Tr20	
	NeonateHead		on	on	100		100None/None	407	69		6,12	92,22	16	standard		4,61	1	192	96	0,6	5	0,8		120	5	Hp38	
InnerEarUHR	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr60	
	InnerEarUHR		on	on	120		120None/None	126	63		8,65	50,34	16	narrow		4,45	1	64	32	0,6	0,4	0,85		40	0,4	Ur77	
InnerEarUHR Sn	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr60	
	InnerEarUHR		on	on	150		150TinFilter/None	473	237		12,92	60,89	16	narrow		10,33	1	64	32	0,6	0,4	0,25		40	0,4	Ur77	
InnerEar	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr60	
	InnerEar		on	on	120		120None/None	126	63		7,97	70,67	16	narrow		3,26	1	192	96	0,6	0,6	0,8		58	0,6	Hr59	
InnerEarSeq	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr60	
	InnEarSeq		on	on	120		120None/None	126	63		8,98	51,72	16	standard		1	1	192	96	0,6	0,6		48	47,7	0,6	Hr59	
Sinus	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Sinus		on	on	100		100None/None	48	24		1,87	16,94	16	narrow		3,3	1	192	96	0,6	3	0,8		60	3	Hr59	
																					3			60	3	Hr40	
																					1			60	0,7	Hr59	
																					1			60	0,7	Hr40	
Orbit	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Orbit		on	on	100		100None/None	120	60		4,67	41,43	16	narrow		3,26	1	192	96	0,6	3	0,8		58	3	Hr59	
																					3			58	3	Hr40	
																					0,6			58	0,4	Hr59	
NeckRoutine	Topogram				150		150TinFilter/None		20		0,01	0,35	32	narrow		2,73		6	6	0,6				256		Tr20	
	Neck		on	on	100		100None/None	163	82		2,6	52,1	32	narrow		5,7	1	192	96	0,6	5	0,8		170	5	Br40	
																					1,5			170	1	Br59	
NeckThorAbd	Topogram				100		100TinFilter/None		30				32	narrow		10,41		6	6	0,6				1024		Tr20	
	Neck		on	on	100		100None/None	163	82		2,6	44,31	32	narrow		5,04	1	192	96	0,6	5	0,8		140	5	Br40	
																					1,5			140	1	Br59	
		Pause																									
	ThorAbd		on	on	120		120None/None	147	74		3,95	257,67	32	narrow		4,72	0,5	192	96	0,6	3	1,4		600	3	Br40	
																					3				3	BI57	
																					3				3	Br40	
																					0,6				0,4	Br40	
ShoulderRoutine	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6				256		Tr20	
	Shoulder		on	on	100		100None/None	178	89		2,82	45,28	32	narrow		4,83	1	192	96	0,6	5	0,8		130	5	Br59	
																					0,75			130	0,5	Br40	
																					0,75			130	0,5	Br59	
TurboFlash_Thorax_HR	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6				512		Tr20	
	FI_Thorax		on	on	100	100	100None/100None	64	32		1,01	20,89	32	narrow	narrow	0,84	0,25	96	96	0,6	3	1,9		150	3	Br40	
																					3			150	3	BI57	
																					0,6			150	0,4	Br40	
																					0,6			150	0,4	BI57	
TurboFlash_Thorax_HR_Sn	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6				512		Tr20	
	FI_Thorax		on	on	100	100	100TinFilter/100TinFilter	408	204		0,71	14,54	32	narrow	narrow	0,84	0,25	96	96	0,6	3	1,9		150	3	Br40	
																					3			150	3	BI57	
																					0,6			150	0,4	Br40	
																					0,6			150	0,4	BI57	
Thorax_XCARE	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6				512		Tr20	
	Thorax		on	on	100		100None/None	64	32		1,01	17,23	32	narrow		1,59	0,25	192	96	0,6	3	0,6		150	3	Br40	
																					3			150	3	BI57	
																					0,6			150	0,4	Br40	

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Child																											
LungLowDose_HR_Sn	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	LungLowDose		on	on	100		100TinFilter/None	101	51		0,18	3,56	32	narrow		0,97	0,25	96	96	0,6	3	1,4			150	3	Br40
																					3				150	3	BI57
																					0,75				150	0,5	Br40
																					0,75				150	0,5	Br64
TurboFlash_Abdomen	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	FI_Abdomen		on	on	100	100	100None/100None	208	104		3,29	84,35	32	narrow	narrow	0,96	0,25	192	96	0,6	5	1,9			200	5	Br40
AbdomenRoutine	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	Abdomen		on	on	100		100None/None	207	104		3,29	83,03	32	narrow		2,24	0,5	192	96	0,6	5	1,4			200	5	Br40
PelvisRoutine	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr60
	Pelvis		on	on	100		100None/None	207	104		3,29	59,49	32	narrow		2,63	0,5	192	96	0,6	5	0,8			150	5	Br40
Spine_Vol	Topogram				150		150TinFilter/None		20		0,01	0,35	32	narrow		2,73		6	6	0,6					256		Tr60
	Spine		on	on	120		120None/None	231	116		6,19	118,08	32	narrow		5,48	1	192	96	0,6	2	0,8			160	2	Br40
																					0,6					0,4	Br40
																					0,6					0,4	Br59
SpineRoutine	Topogram				150		150TinFilter/None		20		0,01	0,35	32	narrow		2,73		6	6	0,6					256		Tr60
	Spine		on	on	120		120None/None	231	116		6,19	118,08	32	narrow		5,48	1	192	96	0,6	2	0,8			160	2	Br40
																					0,6				160	0,4	Br59
ExtrUHR	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	ExtrUHR		on	off	120		120None/None	42	42		2,43	19,03	32	narrow		5,68	1	64	32	0,6	0,4	0,85			60	0,4	Ur77
ExtrUHR_Sn	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	ExtrUHR		on	off	150		150TinFilter/None	110	42		1,11	8,09	32	narrow		7,71	1	64	32	0,6	0,4	0,55			60	0,4	Ur77
Extremity	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	Extremity		on	off	100		100None/None	35	35		1,11	14,48	32	narrow		4,17	1	192	96	0,6	3	0,8			100	3	Br40
																					1				100	0,7	Br40
ExtrUHR	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	ExtrUHR		on	off	120		120None/None	42	42		2,43	19,03	32	narrow		5,68	1	64	32	0,6	0,4	0,85			60	0,4	Ur77
ExtrUHR_Sn	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	ExtrUHR		on	off	150		150TinFilter/None	110	42		1,11	8,09	32	narrow		7,71	1	64	32	0,6	0,4	0,55			60	0,4	Ur77
Extremity	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	Extremity		on	off	100		100None/None	35	35		1,11	14,48	32	narrow		4,17	1	192	96	0,6	3	0,8			100	3	Br40
																					1				100	0,7	Br40
HeadAngio_XCARE	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6					256		Tr20
	HeadAngio		on	on	80		80None/None	155	78		3,21	24,79	16	narrow		1,15	0,25	128	64	0,6	2	0,6			60	2	Hv36
																					0,6				60	0,4	Hv36
HeadAngio	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6					256		Tr20
	HeadAngio		on	on	80		80None/None	155	78		3,21	31,39	16	narrow		1,65	0,5	128	64	0,6	2	1,2			60	2	Hv36
																					0,6				60	0,4	Hv36
CarotidAngio	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	CarotidAngio		on	on	80		80None/None	156	78		1,22	25,4	32	narrow		2,58	0,5	128	64	0,6	5	1,4			170	5	Bv36
																					0,6				170	0,4	Bv40
TurboFlash_BodyAngio	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	FI_BodyAngio		on	on	80	80	80None/80None	156	78		1,14	29,32	32	narrow	narrow	0,96	0,25	192	96	0,6	3	1,9			200	3	Bv36
																					0,6				200	0,4	Bv40
BodyAngio	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	BodyAngio		on	on	80		80None/None	156	78		1,14	28,86	32	narrow		2,24	0,5	192	96	0,6	3	1,4			200	3	Bv36
																					0,6				200	0,4	Bv40
TurboFlash_Cardiac	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	TestBolus	Contrast	off	off	80		80None/None	37	19		3,75	1,87	32	narrow		0,25	0,25	1	1	5	5		0		0		Bv40
	FI_CorCTA		on	on	80	80	80None/80None	444	222		1,02	11,36	32	narrow	narrow	0,59	0,25	192	96	0,6	3	3,2			70	3	Bv40
																					0,6				70	0,3	Bv40
DS_CoroCTA_AdaptSeq	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	DS_CorAdSeq		on	on	80	80	80None/80None	444	222		1,92	20,1	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48		102,6	3	Bv40
																					0,6				102,6	0,3	Bv40
DS_CoronaryCTA	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	DS_CorCTA		on	on	80	80	80None/80None	444	222		2,92	27,79	32	narrow	narrow	2,3	0,25	192	96	0,6	3	0,17			70	3	Bv40
																					0,6				70	0,3	Bv40
																					1,5				70	1	Qr36
DS_CoroCTA_AdaptSeqArrhythmicHR	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	DS_CorAdSeq		on	on	80	80	80None/80None	444	222		4,69	45,05	32	narrow	narrow	0,18	0,25	192	96	0,6	3		48		102,6	3	Bv40
																					0,6				102,6	0,3	Bv40
DS_CoroCTAArrhythmicHR	Topogram				100		100TinFilter/None		30				32	narrow		2,73		6	6	0,6					256		Tr20
	DS_CorCTA		on	on	80	80	80None/80None	444	222		3,83	36,46	32	narrow	narrow	2,3	0,25	192	96	0,6	3	0,17			70	3	Bv40
																					0,6				70	0,3	Bv40
RDE_Thorax_Analysis	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	DE_Thorax		on	off	70	150	70None/150TinFilter	95	48	18	0,87	15,42	32	narrow	narrow	1,69	0,25	192	96	0,6	5	0,55			150	5	Br40
																									150	5	BI64
																					1,5				150	1	Qr40
RDE_Abdomen_VNC	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6					512		Tr20
	DE_Abdomen		on	off	70	150	70None/150TinFilter	170	85	28	1,45	33,3															

Scan Protocol Name	Range Name	Memo	CARE Dose[4D]	CARE kV	Tube Voltage A [kV]	Tube Voltage B [kV]	kV Tin combi	Qual. ref. Tube Current [Qual. ref. mAs]	Tube Current A [mAs]	Tube Current B [mAs]	CTDIvol [mGy]	DLP total [mGy cm]	CTDI Phantom Size [cm]	Shaped Filter A	Shaped Filter B	Scan time (s)	Rot. time [s]	Acquisition	No. Slice Collimation	Coll. Width [mm]	Slice Width [mm]	Pitch Factor	Feed/ Scan [mm]	Length (mm)	Increment (mm)	Kernel	
Child																											
RDE_Spine_Metal_Vol	Topogram				150		150TinFilter/None		20		0,01	0,35	32	narrow		2,73		6	6	0,6				256		Tr20	
	DE_Spine		on	off	100	150	100None/150TinFilter	66	33	66	2,65	52,42	32	narrow	narrow	2,99	0,5	192	96	0,6	3	0,7		160	3	Br40	
																					1,5			160	1	Qr40	
																					1,5			160	1	Br40	
																					2				2	Br40	
																					2				2	Br59	
RDE_Head_BrainHem	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	DE_Head		on	off	70	150	70None/150TinFilter	340	170	74	9,42	140,29	16	standard	standard	3,24	0,5	128	64	0,6	4	0,7		120	4	Hr36	
																					1			120	0,7	Qr40	
RT_Head	Topogram				150		150TinFilter/None		20		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	Head		on	semi	120		120None/None	287	144		18,24	250,54	16	narrow		5,8	1	48	48	1,2	3	0,55		120	3	Hr40	
RT_Abdomen	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6				512		Tr20	
	Abdomen		on	semi	120		120None/None	175	88		4,7	103,49	32	narrow		3,9	0,5	48	48	1,2	3	0,6		200	3	Br40	
TurboFlash_NeonateBody	Topogram				100		100TinFilter/None		20				32	narrow		2,73		6	6	0,6				256		Tr20	
	FI_Neonate		on	on	70	70	70None/70None	448	76		0,68	14,09	32	narrow	narrow	0,84	0,25	192	96	0,6	5	1,9		150	5	Br40	
NeonateBody	Topogram				100		100TinFilter/None		20				32	narrow		2,73		6	6	0,6				256		Tr20	
	NeonateBody		on	on	70		70None/None	448	76		0,68	13,81	32	narrow		1,93	0,5	192	96	0,6	5	1,4		150	5	Br40	
ThoraxUHR_Sn	Topogram				100		100TinFilter/None		30				32	narrow		5,29		6	6	0,6				512		Tr20	
	ThoraxUHR		on	on	150		150TinFilter/None	91	46		1,13	19,32	32	narrow		4,91	0,5	32	32	0,6	3	1		150	3	UI64	
																					3			150	3	Ub36	
																					0,6			150	0,4	UI57	
																					0,6			150	0,4	Ub36	
Sn_Topos		Neonate Sn10																									
	TopoHead				150		150TinFilter/None		55		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	TopoHeadNeck				150		150TinFilter/None		55		0,03	0,72	16	narrow		2,73		6	6	0,6				256		Tr20	
	TopoNeck				150		150TinFilter/None		55		0,01	0,35	32	narrow		2,73		6	6	0,6				256		Tr20	
	TopoNeckThor				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoShoulder				100		100TinFilter/None		82				32	narrow		2,73		6	6	0,6				256		Tr20	
	TopoThorax				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoThorUpAbd				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoTorAbd				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoAbdomen				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoPelvis				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoLegs				100		100TinFilter/None		159				32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoSpine				150		150TinFilter/None		106		0,01	0,69	32	narrow		5,29		6	6	0,6				512		Tr20	
	TopoRunoff				100		100TinFilter/None		312				32	narrow		10,41		6	6	0,6				1024		Tr20	
	TopoWholeBody				100		100TinFilter/None		455				32	narrow		15,17		6	6	0,6				1500		Tr20	

Caution: Federal law restricts this device to sale by or on the order of a physician, dentist, or veterinarian (21 CFR 801.109(b)(1)).

The original language of this document is English.

Siemens Healthineers Headquarters	Legal Manufacturer
Siemens Healthcare GmbH	Siemens Healthcare GmbH
Henkestr. 127	Henkestr. 127
91052 Erlangen	91052 Erlangen
Germany	Germany
Phone: +49 9131 84-0	
siemens-healthineers.com	